

# **NARCH Guidance for water management**

### On a daily basis you should:

- Test the pool/treadmill water 2-3 times per day (more if re-tests are required following adding any chemicals) if the pool/treadmill is operational. The following should be recorded when testing:
  - Free chlorine (bromine)
  - Total chlorine (bromine)
  - o pH
  - o Temperature
  - Who did the test and at what time
  - Any action taken
  - o Re-test results if chemicals added

**Note:** Action taken may include the addition of any chemicals, the quantity of chemicals added, cleaning of pump baskets/skimmer baskets, backwashing etc

#### On a weekly basis you should:

- Backwash or clean filters (more often in pools/water treadmills with a heavy load or low water volume);
- Calculate water balance is water corrosive or scale forming? temperature, hardness, alkalinity, total dissolved solids
- Review test results/patterns and adjust procedures if necessary to maintain correct disinfection levels

## On a monthly basis you should:

- Send a sample of both pool and water treadmill water for microbiological water testing. The testing should include: heterotrophic colony count; coliform bacteria, and pseudomona aeruginosa
- If applicable send a sample for legionnaires disease for more information follow this link: http://www.hse.gov.uk/healthservices/legionella.htm

## On an annual basis you should:

- Service and check all electrical and mechanical equipment
- Calibrate water test equipment
- Change the filter agent

#### **Recommended levels:**

Free Chlorine 2ppm-3ppm – shock to over 6ppm on a weekly basis. Do not treat dogs if under 1ppm

pH: 7.2-7.6

Bromine 4ppm-6ppm

Temp: The safe and comfortable temperature range is 28.0 - 31 degrees celcius. However it is strongly advised to stay at 29 to 30 degrees celcius to obtain maximum therapeutic benefits for

Updated: May 2021



patients. If the pool building is particularly warm or in very hot weather conditions it may be advisable to keep water temperature between 28 and 29 degree Celsius.

## Approved equipment for testing:

- Comparator
- Photometer
- Pool and Spa test strips

Updated: May 2021